TRIM DISC

TRAININGSANLEITUNG // TRAINING INSTRUCTIONS
CONSIGNES D'ENTRAÎNEMENT // GUIDA ALL'ALLENAMENTO

NEW: Trainingsplan, Trainingsvideos & Workouts
www.mft-bodyteamwork.com/trainingsprogramme
(ersetzt DVD/ replaces DVD)

**MFT**[®]
Bodyteamwork

MFT TrimDisc


Übungen

Exercises

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
- Take your time and find a quiet place.
- Train barefoot or with indoor training shoes.
- Find a fixation point on the wall at eye level and look straight ahead. It is best to pay attention to breathing (e.g.: 10 breaths = 10 repetitions = 1 series).
- All exercises should be performed in a standing position with an upright upper body and slightly bent knees. The bobbing and turning movements are supported by a slight movement of the pelvis.
- After a successful exercise: Try the same exercise once more with your eyes closed.
- Choose the exercises and number of repetitions according to your performance level.
- You get the best training effect with approx. 10 – 15 minutes every day.
- If you have health problems or sports injuries, consult your doctor or therapist before the exercises.

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
Move your centre of gravity right and left 10 times, afterwards hold your balance for 10 seconds. Exercise with stretched out knees, lift and lower the hips in turns. The exercise is performed correctly when no "clack!" sound is heard when you touch the disc.

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
Move your centre of gravity forward and backward 10 times. Keep your knees slightly bent! afterwards hold your balance for 10 seconds. The exercise is performed correctly when no "clack!" sound is heard when you touch the disc.

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Do 10 dynamic clockwise turns and 10 anticlockwise by shifting your weight on the right leg and your left leg whilst holding your balance for 10 seconds.

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
Do 10 alternate rhythmic twists without tilting the stand disc. afterwards hold your balance for 10 seconds.

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Lift your hips slowly, hold the disc balanced, then lower again. Repeat 10 times. The exercise is performed correctly when the disc is balanced. Alternative: Tilt the disc forward and backward.

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Move the upper part of the body forward and back 10 times by bending and stretching the arms. Repeat 10 times. Alternative: Tilt the disc whilst keeping it under control.

quiet spot.
 r training shoes.
 wall at eye-level and look
 djust the exercises to your
 ginning (e.g. 10 breaths = 10
 uted in a standing position with
 upright and slightly bent knees;
 ovements are supported by a
 is.
 lose your eyes, breathe deeply
 ce more with closed eyes (to
 .
 mber of repetitions in keeping
 .
 ct when doing the MFT training
 ry day.
 or have suffered sports injuries,
 pist before starting with the

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ght and left 10 times.
 for 10 seconds. Alternative:
 nees while lifting and lowering
 is correct if you do not hear a
 ase plate.

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rward and backward 10 times.
 Afterwards hold your balance
 correct if you do not hear a
 ase plate.

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a-rounds on the stand disc and
 our weight from between your
 le rotating your hips. Afterwards
 nds.

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sts to the left and to the right,
 at the same time. Afterwards
 nds.

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osition, and slowly lower
 xercise is correct if you keep the
 the stand disc forward and

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ody forward by bending your
 stretching them out again.
 Tilt the disc to the left and right
 ol.

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Nach 2001/95/EG
www.mft-bodyteamwork.com/CE

MFT Trainings Devices Instruction Manual

Intended Usage: MFT training devices can be used in athletic and fitness training to improve ones coordination and proprioception and in Physical Therapy for injury prevention and rehabilitation of ones ankles, knees, hips and spine.

Requirements: Constant supervision and instruction from a doctor, certified trainer or therapist is required when using MFT training devices for Physical Therapy.

Age: from 4 to 100 years, according to the level of fitness.

Weight: max. 120 kg body weight (aprox: 270 lbs). Fun Disc 80 kg (aprox: 180 lbs)

Skills: no skills necessary

Safety Instructions: One must execute all workouts on a compact, non slippery, dry and level surface must be clean and dry during all training. One can train with sneakers or barefoot. Do not wear socks, stockings or dress shoes. If necessary, one can use a chair, table, bar, or spotting person for assistance during fitness training. One can reduce the risk of injuries when falling by keeping the workout area clear of debris and if needed, place a soft surface around the training area. Children must only exercise under supervision of an adult.

Additional Instruction: MFT devices must not be used for jumping, but only for smooth and easy movements.

Regulations/Maintenance: Before usage check the device for damages and make sure that all adjustable pieces are fixed.

According to 2001/95/EG
www.mft-bodyteamwork.com/CE

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